


medartis®

PRECISION IN FIXATION

PRODUCT INFORMATION

MODUS® 2
Intermaxillary
Fixation System IMF



MODUS®

Contents

3	MODUS 2 IMF
4	IMF SpeedTip Screws
4	Features
5	IMF Screws with and without Plateau
6	Titanium Arch Bar
7	Technology
9	Storage
10	Ordering Information
12	Bibliography


For further information on the MODUS product line, please visit www.medartis.com.

MODUS 2 IMF

MODUS 2 IMF is indicated for temporary perioperative fixation and/or stabilization of occlusion.

Areas of Use


- Immediate measure to immobilize mandibular fractures or for conservative treatment
- Intraoperative fixation of occlusion
- Postoperative fixation to immobilize the mandible



IMF SpeedTip Screws

Features

- Enlarged head diameter to secure the elastomeric ligatures
- Self-drilling SpeedTip thread design for direct screw insertion without predrilling*
- Reinforced screw core for increased shear strength
- Through bores in which to insert the ligature wires are placed parallel to the markings on the screw head
- Flat head design for soft tissue protection




*With very hard bone it may still be necessary to predrill

IMF Screws with and without Plateau


IMF Screws without Plateau (M2-5248.XX)

The screw head lies directly on the bone, reducing implant protrusion.



IMF Screw with Plateau (M2-5249.xx)


The plateau is over the gingival mucosa and holds wires and elastics away from the soft tissue.



Titanium Arch Bar

Features

- Precontoured pins no longer need to be bent*
- Wavelike arch bar prevents slippage of the ligature wires
- Medial pin serves as an alignment and fixation aid when fitting the arch bar



*If individual pins touch the mucosa, the pins concerned have to be bent away slightly to prevent soft tissue irritation.

Technology

HexaDrive


- A single screwdriver for all screws thanks to the uniform HexaDrive screw head design
 - HexaDrive interface with self-holding properties between screw and screwdriver
 - Increased torque transmission
 - Simplified screw pick-up due to patented self-holding technology



Blue contact surface for torque transmission
Red contact surface for screw retention


SpeedTip

- Functionally unique cutting with immediate bite
- Immediate cutting of the bone with only slight axial pressure
- The triangular tip design permits simultaneous drilling, tapping and compression of the bone tissue during insertion for increased pull-out stability ⁽¹⁻²⁾
- Reduced insertion torque thanks to the polygonal tip and tapered shaft



Storage

- Compact, clearly arranged system
- Easy to use
- User-friendly storage of implants and instruments
- Validated cleaning and sterilization of the implants in the container



Ordering Information

2.0 IMF SpeedTip Screws, HexaDrive 6

Material: Titanium alloy (ASTM F136)



Length	Art. No.	STERILE	Description	Pieces/Pkg
8 mm	M2-5248.08	M2-5248.08S	without plateau	2
11 mm	M2-5248.11	M2-5248.11S	without plateau	2
14 mm	M2-5248.14	M2-5248.14S	without plateau	2

2.0 IMF Screws with Plateau, HexaDrive 6


Material: Titanium alloy (ASTM F136)



Length	Art. No.	STERILE	Description	Pieces/Pkg
8 mm	M2-5249.08	M2-5249.08S	with plateau	2
11 mm	M2-5249.11	M2-5249.11S	with plateau	2

Titanium Arch Bar

Material: Titanium alloy (ASTM F67)
Plate thickness: 0.7 mm



Art. No.	STERILE	Length	Pieces/Pkg
M-4450	M-4450S	129 mm	2

Twist Drills Ø 1.5 mm (Core Hole 2.0 Screws)



M2-3459



M2-3469

Art. No.	STERILE	Description	Stop	Length	Shaft End	Pieces/Pkg
M2-3459	M2-3459S	for drill guide M2-2198	25 mm	99 mm	Dental	1
M2-3469	M2-3469S	for drill guide M2-2198	25 mm	112 mm	Stryker J-Latch	1

Screwdriver Handle



1:2

M2-2001

Art. No.	Description	Length	Pieces/Pkg
M2-2001	Type 2 (hand-driven small, AO coupling)	121 mm	1

Screwdriver Blade



Art. No.	Interface	Description	Length	Pieces/Pkg
M2-2005	⊙ HD6	self-holding	95 mm	1

Bibliography

Self-drilling Screws

1. Heidemann, W.; Terheyden, H.; Gerlach, K. L.
Analysis of the osseous / metal interface of drill free screws and self-tapping screws
Journal of Cranio-Maxillofacial Surgery
(2001) 29, 69–74
2. Heidemann, W.; Terheyden, H.; Gerlach, K. L.
In-vivo-Untersuchungen zum SchraubenKnochen-Kontakt von Drill-Free-Schrauben und herkömmlichen selbstschneidenden Schrauben
Mund Kiefer GesichtsChir 5 2001: 17–21

MANDIBLE2-01020001_v0 / © 2020-07, Medartis AG, Switzerland. All technical data subject to alteration.

MANUFACTURER & HEADQUARTERS

Medartis AG | Hochbergerstrasse 60E | 4057 Basel/Switzerland
P +41 61 633 34 34 | F +41 61 633 34 00 | www.medartis.com

SUBSIDIARIES

Australia | Austria | Brazil | France | Germany | Japan | Mexico | New Zealand | Poland | UK | USA

For detailed information regarding our subsidiaries and distributors, please visit www.medartis.com



Disclaimer: This information is intended to demonstrate the Medartis portfolio of medical devices. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. For further questions, please contact your Medartis representative (www.medartis.com). This information contains CE-marked products.
For US only: Federal law restricts this device to sale by or on the order of a physician